

Sterilization Solutions



Euro Clinical Service



DuPont™
Tyvek®



Linea
STERIPLASMA

Pouches and reels in Tyvek®*

Reels and pouches are made of a transparent plastic film in polyester and polyethylene named PET12/PE50, and of a TYVEK® side of the company DuPont™, with colour changing indicator of process certified ISO 11140-1.

The STERIPLASMA LINE is suitable to Plasma sterilization with Hydrogen peroxide VH2O2, to Ethylene oxide EO sterilization and to radiation sterilization IRRAD.

Raw materials used to produce the STERIPLASMA LINE grant a high resistance degree to laceration and a high barrier to bacteria.

The complete STERIPLASMA LINE is produced according to a such high quality standard level to grant a perfect peelability and avoid the release of fibers or particles on the surfaces.

All cartons are endowed with use instructions in various languages.

Tyvek® is a registered mark of the company Du Pont™.

Buste e rotoli in Tyvek®*

I rotoli e le buste sono composti da un film plastico trasparente in poliestere e polietilene denominato PET12/PE50 e da un lato in TYVEK® della società DuPont™, con indicatore di processo a cambiamento di colore certificato ISO 11140-1.

La LINEA STERIPLASMA è idonea alla sterilizzazione al plasma con perossido di idrogeno VH2O2, alla sterilizzazione ad ossido di etilene EO ed irraggiamento IRRAD.

Le materie prime utilizzate per la produzione della LINEA STERIPLASMA garantiscono un'elevata resistenza alla lacerazione ed un'elevata barriera ai batteri.

Tutta la LINEA STERIPLASMA è prodotta seguendo un elevato standard qualitativo tale da garantire una perfetta pelabilità ed il non rilascio di fibre e particelle sulle superfici.

Tutti i cartoni contengono istruzioni d'uso in varie lingue.

Tyvek® è un marchio registrato Du Pont™.



Linea STERIPLASMA



DuPont™
Tyvek®

FLAT REELS ROTTOLI PIATTI TYVEK®

Ref.	Dimensions
TVRP7570	75mmx70mt
TVRP1070	100mmx70mt
TVRP1570	150mmx70mt
TVRP2070	200mmx70mt
TVRP2570	250mmx70mt
TVRP3070	300mmx70mt
TVRP3570	350mmx70mt
TVRP4270	420mmx70mt
TVRP5070	500mmx70mt

other sizes upon request
altri formati a richiesta

FLAT POUCHES BUSTE PIATTE TYVEK®

Ref.	Dimensions
TVBP7520	75mmx200mm
TVBP1026	100mmx260mm
TVBP1532	150mmx320mm
TVBP2040	200mmx400mm
TVBP2050	200mmx500mm
TVBP2550	250mmx500mm
TVBP3065	300mmx650mm

other sizes upon request
altri formati a richiesta

Self-Sealing Pouches
upon request

Buste auto saldanti
a richiesta